

## 4 Point Bend Fixture for Glass Testing | CP117856

### 4 Point Bend Fixture for Glass Testing

The 4 point bend fixture for glass testing allows for the testing of glass samples up to 30cm x 30cm. It features 3 articulating and rolling anvils which utilize needle roller bearings for repeatable rolling from test to test even under load. The maximum capacity of this fixture is 200N.

### Specifications

Maximum Load	200 N
Base Connection	Type Dm
Load Cell Connection	Type Om
Machine Compatibility	Dual Column Machines
Maximum Upper Span	200mm
Maximum Lower Span	300mm

