

2kN High Temp Pneumatic Side Action Grips | CP122505

2kN High Temp Pneumatic Grips

2kN pneumatic side action grips designed for use in non-ambient temperature ranges

- Based on the 2712-041 pneumatic grips. Compatible with accessories and jaw faces
- Stainless steel construction for use at high and low temperatures
- Integrated pullrods
- Grips can be positioned at 0, 30, and 90 degrees for use with AVE2/SVE2
- Requires pneumatic control kit

Specifications

Max Load	2 kN
Max. Specimen Thickness	13 mm (0.5 in)
Connection Type	6 mm clevis (Type 0m)
Temperature Range	-73°C to 316°C (-100°F to 600°F)
Grip Separation with 3119-609 chamber on 596X E2 Frame	0 mm to 390 mm



www.instron.com